



TIMEX SENSOR

THE EVOLUTION OF YOUR LOOK

It's time to get out of the analog period and scan your clinic for clearer, more accurate and safer diagnosis.



Size 1
Active area: 20 X 30 mm

Size 2
Active area: 26 X 36 mm

Timex Sensor offers powerful imaging and connects your office to the future.

Using it you optimize intraoral radiographic capture processes and efficiently manage your patients images.

Timex Sensor has intuitive software that is compatible with periapical x-ray equipment and can be installed on as many computers as you wish

Forget radiographic films, developers and fixers. Save time and contribute to the preservation of the environment!

BENEFITS

- Facilitates access among professionals
- Manipulate and save images
- Environmentally-friendly
- Provides accuracy in diagnosis
- Reduce film and processing costs

Features

Replaceable Cable

Timex Sensor uses an easy-to-replace 3 m USB cable.

Shock resistant

The internal protective silicon cover makes Timex shock resistant.

Unlimited licenses

There are no license limits and you can install Timex Sensor Software on any computer* in your clinic, save the images and access them from anywhere.

The perfect Picture sharpness

Capture high definition images in real time. With 25 line pairs (theoretical) and advanced image processing features, you have clear images for the best diagnosis.

Organized images

Register your patients, capture images instantly and organize all images in seconds!

Share images and diagnostics

Save images in various formats like TIFF, BITMAP, PNG, DICOM, JPEG. Generate PDF reports and instantly share them with your clientes and colleagues who will continue the treatment.

Full feature

Wide range of tools such as rotation, zoom, brightness, contrast and gamma adjustment, positive/negative, color adjustment, embossing, text insertion, arrows, circles, area selection and cropping, angle and segment measurement, bar for reports/diagnosis, odontogram and comparison of full screen images.

Save time

While it takes you about 5 minutes to capture and reveal an image through the analog system, Timex Sensor does it in about 5 seconds!

Compatibility

Timex Sensor Intraoral Sensor is compatible with existing intraoral X-ray equipment.

Save Money

Stop spending on movies developers, fixers and chemical disposal rate. The investment to have your Timex Sensor is less than amount spent monthly on these inputs.

Impress your customers

Your clients will be impressed by the technology and richness of detail before and after the procedure.

Ecologically correct

Stop discarding chemicals as developers and fixers.

Technical Information

Operational systems: Windows 7 Home, Start, Pro (x86 and x64); Windows 8 Pro (x86 and x64); Windows 8.1 Pro (x86 and x64); Windows 10 Home e Pro (x86 and x64) e CPU: Intel Core i3 3.0 GHz or higher. RAM: 4GB DDR2 SDRAM- HDD 500GB or higher- USB Connector: 2x USB 2.0 – Monitor 1024 X 768.



SAEVO GET YOURS!